



Università degli Studi di Milano
Department of Pharmaceutical Sciences



Speech commands

Alessandro Pedretti

LabSim - Speech commands

Command	Event
Turn on (the) bunsen Turn on (the) burner	If you selected the lighter, the burner ignites.
Turn on (the) centrifuge	The centrifuge turns on.
Turn on (the) light	The lights comes on :-D
Turn on (the) TV	The TV turns on.
Smell	Smells the selected glassware. If the glassware is not selected, smell the air.
Smell (the) stirrer	Smells the stirrer.
Smell (the) beaker	Smells the beaker.
Smell (the) dropper	Smells the dropper.
Smell (the) watch glass	Smells the watch glass.
Smell (the) curcible	Smells the curcible.
Smell (the) test tube	Smells the selected test tube. If no test tube is selected, smell the tube number 1.
Smell (the) test tube (number) 1	Smells the test tube number 1.
Smell (the) test tube (number) 2	Smells the test tube number 2.
Smell (the) test tube (number) 3	Smells the test tube number 3.
Smell (the) test tube (number) 4	Smells the test tube number 4.
(Run) (the) centrifuge for one minute	Runs the centrifuge for one, two, three or four minutes.
(Run) (the) centrifuge for two minutes	
(Run) (the) centrifuge for three minutes	
(Run) (the) centrifuge (for) four minutes	
Open (the) centrifuge	Opens the centrifuge lid if it is previously unlocked.
Open (the) (bunsen) tap	Opens the bunsen tap.
Close (the) centrifuge	Closes the centrifuge lid.
Close (the) (menu)	Closes the menu currently opened.
Close (the) (bunsen) tap	Closes the bunsen tap and turn off the flame if lit.
Oxidizing flame	Sets the bunsen with the oxidizing flame.
Reducing flame	Sets the bunsen with the reducing flame.
Turn (the) ring	Turns the bunsen ring changing the flame.
About	Shows the information window.
Leave	Leaves the glassware.
Leave (the) lighter	Leaves the lighter.
Leave (the) stirrer	Leaves the stirrer.
Leave (the) beaker	Leaves the beaker.
Leave (the) dropper	Leaves the dropper.
Leave (the) watch glass	Leaves the watch glass.
Leave (the) curcible	Leaves the curcible.
Leave (the) test tube	Leaves the test tube.
Leave (the) test tube (number) 1	Leaves the test tube number 1.
Leave (the) test tube (number) 2	Leaves the test tube number 2.
Leave (the) test tube (number) 3	Leaves the test tube number 3.
Leave (the) test tube (number) 4	Leaves the test tube number 4.
Menu Main menu	Shows the main menu.

LabSim - Speech commands

Analyte menu	Shows the analyte menu.
Settings menu	Shows the settings menu.
Settings	
Put (it) in the water bath	Puts the selected glassware in the water bath.
Put (ii) in the centrifuge	Puts the selected test tube in the centrifuge.
I don't understand I didn't understand What	:-D
Ok Google	:-D
Take (the) lighter	Takes the lighter.
Take (the) stirrer	Takes the stirrer.
Take (the) beaker	Takes the beaker.
Take (the) watch glass	Takes the watch glass.
Take (the) curcible	Takes the curcible.
Take (the) test tube	Takes the test tube.
Take <RACTIVE_NAME>	The reactive is selected and moved to foreground. Synonyms of the reactive are allowed.
Take analyte	The analyte is selected and moved to foreground.
Unlock (the) centrifuge	Removes the security lock from the centrifuge lid.
Turn off (the) bunsen	The burner turns off.
Turn off (the) centrifuge	The centrifuge turns off.
Turn off (the) light	The light switches off :-D
Turn off (the) TV	The TV turns off.
Use (the) lighter with (the) bunsen Use (the) lighter with (the) burner	Uses the lighter with the bunsen to ignite it.
Use (the) lighter with (the) curcible	Uses the lighter with the curcible to ignite it.
Use (the) stirrer	Uses the stirrer.
Use (the) dropper Use (the) pipette	Adds the liquid reactive with the dropper or uses the generic dropper.
Use (the) (nickel chrome) wire with (the) (test) tube	Dips the nickel-chromium wire in the first tube containing concentrated hydrochloric acid.
Use (the) (nickel chrome) wire with (the) (test) tube (number) 1 Use (the) (nickel chrome) wire with (the) (test) tube (number) 2 Use (the) (nickel chrome) wire with (the) (test) tube (number) 3 Use (the) (nickel chrome) wire with (the) (test) tube (number) 4	Dips the nickel-chromium wire in the specified tube (from 1 to 4) containing concentrated hydrochloric acid.
Use (the) spatula	Uses the spatula.
Use (the) wash bottle	Uses the wash bottle.
Go to (the) workbench	Looks at the workbench.
Go to (the) centrifuge	Looks at the centrifuge.
Go to (the) TV	Looks at the TV.

* The optional words are between parentheses.